

# TMS Targets

## - Overview -



## Disclaimer

All content presented is no medical advice but for informational and educational purposes only. We focus on scientific evidence, not regulatory approvals. The Academy does not recommend or endorse any specific practitioner, equipment, product, medical procedure, or protocol presented in our courses or on our webpages and the Academy does not give direct medical advice. Reliance on any information provided or mentioned during our courses, our recordings or on our websites is always at your own risk. Practitioners should always consult additional sources of information and exercise their best professional judgement before making clinical decisions of any kind.

## TMS Targets

Many of the brain regions discussed as potential **targets for TMS** are **ideally defined anatomically or functionally** in terms of individual structural or functional brain scans.


Such fMRI guided TMS neuronavigation is the most accurate method of targeting a given brain region in a given patient with high accuracy and reliability.

However, the exact added value for clinical outcomes remains debated.

**In our target overviews we discuss simpler and more cost-effective alternatives for targeting certain brain regions, guiding the coil positioning using scalp based landmark approaches.**

Please note that such one size fits all approaches inherently suffer from inter individual variability in the correspondence between the proposed scalp landmark and the assumed underlying brain structures.

## TMS Targets



## Primary motor cortex (M1)

### Coil position/orientation/type:

- **Motor hot spot** (hand area) determination with cortical mapping procedure
- **Coil handle points in the lateral-posterior direction at 45 degrees** relative to the midline
- Standard **figure-eight coil**

### Treatment examples:


#### Left & right M1

Motor stroke (hand recovery)

Neuropathic pain

Fibromyalgia

Parkinson's disease (leg area)



# Dorsolateral prefrontal cortex (DLPFC)


## 5/6/7cm-rule

### Coil position/orientation/type:

- 5/6/7cm anterior to motor hotspot
  - Parallel to midline
- Coil handle points in the lateral-posterior direction at 45 degrees relative to the midline
- Standard **figure-eight coil**

### Treatment examples:

F3	F4
Depression	Depression
Schizophrenia (negative symptoms)	OCD PTSD



## Dorsolateral prefrontal cortex (DLPFC)


### F3/F4 (10/20 system)

**Coil position/orientation/type:**

- **F3/F4** (based on 10/20 system)
- **Coil handle points in the lateral-posterior direction at 45 degrees relative to the midline**
- same as motor cortex stimulation
- Standard **figure-eight coil**

**Treatment examples:**

F3	F4
Depression	Depression
Schizophrenia (negative symptoms)	OCD PTSD




# Dorsolateral prefrontal cortex (DLPFC)

## Beam F3/F4

⚠ Beam F3/F4 is not identical to F3/F4 (10/20 system)

### Coil position/orientation/type:

- Beam F3/F4
- Coil handle points in the lateral-posterior direction at 45 degrees relative to the midline
- same as motor cortex stimulation
- Standard **figure-eight coil**




(F3) Distance along circumference from midline (X):	(F3) Distance from vertex (Y):
<b>6.45 cm</b>	<b>9.46 cm</b> Adjusted*: 9.81 cm

### Treatment examples:

F3	F4
Depression	Depression
Schizophrenia (negative symptoms)	OCD
Addiction	PTSD

**BA9 BA8 BA43 Location System**  
Will Beam & Jeff Borckardt  
Web Interface Developed 6/7/2010

1. Tragus to Tragus (CM):	<input type="text" value="37"/>
2. Nasion to Inion (CM):	<input type="text" value="35"/>
3. Circumference (CM):	<input type="text" value="56,5"/>
<input type="button" value="Calculate"/>	



### Beam F4:

Same coordinates on right hemisphere



## Dorsomedial prefrontal cortex (DMPFC)

### Coil position/orientation/type:

- Move 25% of the nasion-inion distance from nasion to CZ
- Coil handle points backwards along the midline = bilateral DMPFC
- Coil handle points to the right = preferential stimulation of left hemisphere
- Double-cone coil

### Treatment examples:

Bilateral DMPFC	Preference left
OCD	Depression



# Right Orbitofrontal Cortex (OFC)

## AF8 electrode position


### Coil position/orientation/type:

- AF8
  - Move 10% of the nasion-inion distance from the nasion to CZ
  - Move 10% of the head circumference to the right hemisphere
- Coil handle pointing downwards
- Double-cone coil

### Treatment examples:

**Right OFC**

- OCD
- Depression



## Pre-Supplementary Motor Area (pre-SMA)


### Coil position/orientation/type:

- Move 35% of the nasion-inion distance from the nasion to CZ
- Coil handle points to the back
- Standard **figure-eight coil**

### Treatment examples:

**Bilateral pre-SMA**

OCD



# Temporoparietal junction (TPJ) T3/P3 or T4/P4 (10/20 system)


## Coil position/orientation/type:

- Midpoint between T3/P3 or T4/P4 (10/20 system)
- No consensus on preferred coil orientation → **coil handle pointing upwards** seems best
- Standard **figure-eight coil**

## Treatment examples:

**TPJ**

- Tinnitus
- Schizophrenia (positive symptoms)



# Auditory cortex


## T3-T5 / T4-T6 (10/20 system)

### Coil position/orientation/type:

- Midpoint between T3-T5 or T4-T6 (10/20 system)
- Accurate coil positioning may require MRI data
- Coil handle points to the back
- Standard figure-eight coil

**Left auditory cortex**

Tinnitus




**Left auditory cortex**

## Broca's Area

### F7 electrode position

#### Coil position/orientation/type:

- **F7**
- Accurate coil positioning may require the induction of speech arrest or MRI data (BA44/45)
- **Coil handle points to the back**
- Standard **figure-eight coil**



#### Treatment examples:

**Broca's Area**

Alzheimer's disease

# Wernicke's Area

## T7 electrode position


### Coil position/orientation/type:

- T7
- Accurate coil positioning may require MRI data
- **Coil handle points to the back**
- Standard **figure-eight coil**

### Treatment examples:

**Wernicke's Area**

Alzheimer's disease



# Precuneus

## CPz electrode position


### Coil position/orientation/type:

- CPz
- Accurate coil positioning may require MRI data
- **Coil handle points to the back**
- Standard **figure-eight coil**

### Treatment examples:

**Precuneus**

Alzheimer's disease



# Posterior Parietal Cortex (PPC)

## Contralesional P3/P4


### Coil position/orientation/type:

- **Contralesional P3/P4 (10/20 system)**
- **Coil handle points in the lateral-posterior direction at 45 degrees relative to the midline**
- Standard **figure-eight coil**

### Treatment examples:

**PPC**

Hemispatial neglect



# Inferior frontal gyrus (IFG) F7/F8 (10/20 system)


## Coil position/orientation/type:

- F7/F8 (10/20 system)
- Accurate coil positioning may require MRI data (BA45)
- **Coil handle points to the back**
- **Standard figure-eight coil**

## Treatment examples:

**(right) IFG**

Post-stroke aphasia



# Superior temporal gyrus (STG) T3/T4 (10/20 system)


## Coil position/orientation/type:

- T3/T4 (10/20 system)
- Accurate coil positioning may require MRI data
- **Coil handle points to the back**
- Standard **figure-eight coil**

## Treatment examples:

**STG**

Schizophrenia (auditory-verbal hallucinations)



# Inferior parietal lobule (IPL) P3/P4 (10/20 system)


## Coil position/orientation/type:

- P3/P4 (10/20 system)
- Accurate coil positioning may require MRI data
- **Coil handle points in the lateral-posterior direction at 45 degrees relative to the midline**
- Standard **figure-eight coil**

## Treatment examples:

**STG**

Alzheimer's disease



# TMS Targets

## - Overview -